

Maximizing Antibacterial Activity of Polymyxin B: Microbiological component analyses alone and in combination

Zahra Kassamali¹, Paul R. Rhomberg², Amy A. Watters², Randall A. Prince³, Larry H. Danziger⁴, John C. Rotschafer⁵, Ronald N. Jones^{2,6}
University of California, Los Angeles¹, JMI Laboratories², University of Houston³, University of Illinois, Chicago⁴, College of Pharmacy University of Minnesota⁵, Tufts University School of Medicine⁶

Abstract

Objectives: The individual components that comprise polymyxin B vary in their antimicrobial activity. As a product of fermentation, polymyxin B composition varies between manufacturers and lots. One component, polymyxin B1, comprises two-thirds of the drug yet demonstrates the least *in vitro* antimicrobial activity. Our objective was to evaluate which combination of polymyxin B components demonstrates the greatest antibacterial activity.

Methods: Polymyxin B1, B2, Ile-B1, and B3 were obtained from the commercial manufacturer, Toku-E. All combinations of 2 components were evaluated using checkerboard panels. Broth microdilution was conducted according to CLSI M07-A9 with dilutions ranging from 0.06 – 4 µg/ml for each component. 12 strains of lactose- and non-lactose fermenting Gram-negative bacilli were tested in duplicate. Synergy was defined as ≥ 4-fold lowering of both component MICs when combined. Partial synergy was defined as ≥ 4-fold decrease in MIC of one component and a 2-fold decrease in MIC of the 2nd component. Additive was defined as a 2-fold decrease in MICs of both components. Indifference was defined as no change in component MICs. Antagonism was defined as ≥ 4-fold increase in both component MICs.

Results: One combination met criteria for synergy against the Enterobacteriaceae: B3 + Ile-B1. Against *P. aeruginosa* and *A. baumannii*, B3 + Ile-B1 demonstrated partial synergy. No combinations demonstrated antagonism. Combinations employing B3 and/or Ile-B1 demonstrated the greatest rates of partial synergy. The B1/B2 combination was least likely to show any level of synergy (13%), with 42% occurrences of indifference. This contrasts to 0-8% indifference and 42-67% partial or complete synergy for all non-B1 combinations.

Conclusion: The dominant components of polymyxin B products (B1 and B2) are those associated with the lowest probability of improved activity when combined. The minor components demonstrate the greatest antimicrobial activity. Refining the composition of polymyxin B to increase the abundance of more active minor components could enhance antibacterial properties of this pharmaceutical.

Introduction

- Polymyxin B is produced by a fermentative process; composition varies between manufacturers and lots.^{1,3}
- We have demonstrated previously that individual components comprising polymyxin B vary in their antimicrobial activity. The more abundant components demonstrated less *in vitro* antimicrobial activity compared to less abundant components.⁴
- The objective of this investigation was to evaluate antimicrobial activity of polymyxin B components in combination.

Methods

- Polymyxin B1, B2, Ile-B1, and B3 were obtained from the commercial manufacturer, Toku-E (715 W. Orchard Dr. Ste 3 Bellingham, WA 98225 USA).
- Each combination of 2 components was evaluated using checkerboard panels conducted in cation-adjusted Mueller-Hinton broth without 0.002% P-80.
- Dilutions ranged from 0.06 – 4 µg/mL for both drug components
- Twelve strains of Gram-negative clinical pathogens were tested in duplicate.
- Combination activity was defined according to the following categories:

| | |
|------------------------|---|
| Synergy | ≥ 4-fold decrease in each component MIC when combined |
| Partial synergy | ≥ 4-fold decrease in MIC of one component and a 2-fold decrease in MIC of the other component |
| Additive | 2-fold decrease in the MICs of both components |
| Indifference | No change in component MIC values when combined |
| Antagonism | ≥ 4-fold increase in both component MIC values when combined |

Results

Table 1: Activity of polymyxin B components (four) tested in combination

| Combination | Organism Group | Antagonism | Indifference | Additive | P. Synergy | Synergy |
|-------------|-----------------------------------|------------|--------------|----------|------------|---------|
| B1/B2 | Enterobacteriaceae ¹ | 0 | 1 | 9 | 2 | 0 |
| | <i>P. aeruginosa</i> ² | 0 | 6 | 0 | 0 | 0 |
| | <i>A. baumannii</i> ² | 0 | 3 | 2 | 1 | 0 |
| B1/Ile-B1 | Enterobacteriaceae | 0 | 3 | 7 | 2 | 0 |
| | <i>P. aeruginosa</i> | 0 | 0 | 5 | 1 | 0 |
| | <i>A. baumannii</i> | 0 | 0 | 2 | 4 | 0 |
| B1/B3 | Enterobacteriaceae | 0 | 3 | 4 | 5 | 0 |
| | <i>P. aeruginosa</i> | 0 | 4 | 0 | 2 | 0 |
| | <i>A. baumannii</i> | 0 | 0 | 4 | 2 | 0 |
| B2/Ile-B1 | Enterobacteriaceae | 0 | 0 | 6 | 6 | 0 |
| | <i>P. aeruginosa</i> | 0 | 0 | 6 | 0 | 0 |
| | <i>A. baumannii</i> | 0 | 0 | 2 | 4 | 0 |
| B2/B3 | Enterobacteriaceae | 0 | 1 | 5 | 6 | 0 |
| | <i>P. aeruginosa</i> | 0 | 0 | 0 | 6 | 0 |
| | <i>A. baumannii</i> | 0 | 1 | 1 | 4 | 0 |
| B3/Ile-B1 | Enterobacteriaceae | 0 | 1 | 4 | 6 | 1 |
| | <i>P. aeruginosa</i> | 0 | 0 | 5 | 1 | 0 |
| | <i>A. baumannii</i> | 0 | 0 | 1 | 5 | 0 |

¹Enterobacteriaceae: 3 strains *E. coli*, 3 strains *K. pneumoniae*
²Three strains each of *P. aeruginosa* and *A. baumannii*

Figures 1a – 1d: Component activity in combination by organism




Figure 1a




Figure 1b




Figure 1c




Figure 1d

ENT: Enterobacteriaceae, PSA: *Pseudomonas aeruginosa*, ACB: *Acinetobacter baumannii*

Summary

- Synergy was seen with one combination, B3/Ile-B1, against the Enterobacteriaceae
- Synergy of any combinations against either *P. aeruginosa* or *A. baumannii* was not identified
- No instances of antagonism were identified
- Combinations including B1 demonstrated the highest rates of indifference
- Combinations including B3 demonstrated the highest rates of partial synergy

Conclusions

- Our analyses indicate the dominant components of polymyxin B products (B1 and B2) are those associated with the lowest probability of improved antibacterial activity when combined.
- While these results should be confirmed with additional investigations including time-kill studies, preliminary data indicate the formulation of polymyxin B could be refined to maximize antimicrobial activity.

References

- Kassamali Z, Rotschafer JC, Jones RN, Prince RA, Danziger LH. Polymyxins: wisdom does not always come with age. Clin Infect Dis 2013; 57(6):877-83.
- Kassamali Z, Wang W, Jones RN, Rotschafer JC, Danziger LH, Prince RA. Is what you see what you get? Product variation in clinical and reference formulations of polymyxin B. European Congress of Clinical Microbiology and Infectious Diseases, Berlin, Germany, April 2013. Abstract #1119.
- Orwa JA, Govaerts C, Busson R, Roets E, Van Schepdael A, Hoogmartens J. Isolation and structural characterization of polymyxin B components. J Chromatog A 2001;912:369-373.
- Kassamali Z, Rhomberg PR, Jones RN, Danziger LH, Prince RA. Polymyxin B: Can we create a better product? European Congress of Clinical Microbiology and Infectious Diseases, Barcelona, Spain, May 2014. Abstract #P1730.